



## ST Zirconia Block

- Super translucent
- Suitable for full contour crown , bridge and anterior

### Indications



### Physical Characteristics

|                       |                  |
|-----------------------|------------------|
| Sintered density      | 6.07±0.03g/cm³   |
| Bending Strength      | ≥1200MPa(Av.)    |
| Transmittance         | 43%              |
| Hardness              | 1200HV           |
| Sintering temperature | Recommend 1530°C |



## HT Zirconia Block

- High translucent
- Brilliant bending strength and economical Block
- Suitable for coping ,framework

### Indications



### Physical Characteristics

|                       |                  |
|-----------------------|------------------|
| Sintered density      | 6.07±0.03g/cm³   |
| Bending Strength      | ≥1200MPa(Av.)    |
| Transmittance         | 40%              |
| Hardness              | 1200HV           |
| Sintering temperature | Recommend 1530°C |

