

NOVELTIES IDS 2019

NEUHEITEN IDS 2019

MORSA CAD Wax gingiva - Prothetik

Article.Number: DR000732HD830M

For making functional try ins and individual bite rims in the complete denture technique

The process supporting CAD discs are used to transform conventional working steps in the laboratory into complete digital processes.

The discs can also be used for conventional finishing procedures in injection and pressing techniques - the ability to boil it out for wax removal makes it possible



MORSA CAD Wax gingiva

Artikelnummer: DR000732HD830M

MORE EFFICIENCY

WITH MORSA

- Zur Herstellung von Wachs-Einproben und individuellen Wachsbiss Schablonen in der Totalprothetik. Ebenso geeignet für die konventionelle Fertigstellung in der Injektions- und Presstechnik

INDIKATIONEN

- Wachsbiss Schablonen Totalprothetik
- Funktionseinproben

MORSA Gingiva Extra Tough Plattenwachs, innovative Rezeptur 500 g

Artikelnummer: DM000154VA804

Ein hartes Vier-Jahreszeiten-Aufstellwachs mit hoher Erweichungstemperatur. Zur Aufstellung künstlicher Zähne für prothetische Arbeiten, z. B. für partielle und totale Kunststoffprothesen, Modellgussättel, Interimsprothesen und Erweiterungen gleichermaßen geeignet. Die natürlich wirkende gingiva Zahnfleischfärbung harmoniert mit heiß- und kaltpolymerisierenden Prothesenkunststoffen

- zur sicheren Fixierung von künstlichen Zähne
- sehr Hohe Verflüssigungstemperatur bewirkt die Möglichkeit längerer Anproben und die Verwendung von diversen Dupliergelen
- überaus geringe Rückfederung (Flow)
- Höchste Dimensionstreue
- Geringste Erstarrungskontraktion sowie Verzug der Zahnaufstellung nach dem Abkühlen
- saubere Modellierbarkeit + Polierbarkeit
- Natürliche rosa gingiva Farbe korrespondierend zu Prothesenkunststoffen
- Typ 3 Wachs



MORSA Gingiva Extra Tough Sheet wax – new innovative formulation 500 g

Article.Number: DM000154VA804

A tough, four-seasons setting-up wax with high softening temperature. Used for setting up denture teeth in wax for all prosthetic restorations, e.g. equally suitable for partial and full acrylic dentures, chrome cobalt saddles, immediate and temporary dentures and extensions. The natural-looking gingival shading corresponds to hot and cold curing denture acrylics.

- High flow temperature enables longer try-ins and use of duplicating gels
- Minimum spring-back (flow)
- High shape and dimensional accuracy
- Low setting contraction, minimum warpage of the set-up after cooling



MORE EFFICIENCY
WITH MORSA

Morsa Wachswarenfabrik Sallinger GmbH

Nordstraße 3

D-86381 Krumbach

Deutschland

Telefon: + 49 (0) 8282 800 44-0

Fax: + 49 (0) 8282 800 44-33

www.morsa-dental.de

info@morsa.de